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MAINTENANCE HANDBOOK

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## ORIGINAL INSTRUCTIONS

### Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

---

(Applicable in the European Union and other European countries with separate collection system)

This marking shown on the product or its literature, indicates that it should not be disposed with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material sources.

Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

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# CONTENTS

---

1.0	OPERATING CONDITIONS	15
1.1	General information	
1.2	Operating conditions	
2.0	SPECIFICATIONS	15
3.0	INSTALLATION	16
3.1	Introduction	
3.2	Damage check	
3.3	Cleaning	
3.4	Working conditions	
3.5	Levelling	
3.6	Putting into operation	
4.0	CONTROL PANEL	17
5.0	OPERATION	18
5.1	Switching on	
5.2	Switching off	
5.3	Set point	
5.4	Editing the setpoint	
5.5	Work cycle	
6.0	SENSORS	19
7.0	DEFROSTING	19
7.1	Automatic defrosting	
7.2	Manual defrosting	
8.0	CONDENSATION EVAPORATION	19
9.0	MINIMUM AND MAXIMUM TEMPERATURE	19
9.1	Viewing maximum temperature recorded	
9.2	Viewing minimum temperature recorded	
9.3	Erasing recorded values	
10.0	OPERATING TIPS	20
11.0	ALARMS AND WARNING	20
11.1	List of alarms	
12.0	CONNECTION DIAGRAM	21
13.0	MAINTANCE AND CLEANING	21
13.1	Cleaning the appliance	
13.2	Cleaning the condenser	
13.3	Cleaning the gaskets	
14.0	GAS SAFETY DATA SHEET	22
14.1	R290	
14.2	R452A	
15.0	APPENDIX A-B-C-D-E-F	50

# 1.0 OPERATING CONDITIONS

## 1.1 GENERAL INFORMATION

Our appliances are designed and optimised in order to obtain high performance. This appliance must not be used by people with limited physical, mental or sensory abilities or without experience and knowledge of it, unless instructed in its use by those responsible for their safety.

## 1.2 OPERATING CONDITIONS

-  **Do not store pressurised spray bottles bearing the wording "flammable". Risk of explosion!**
-  **Keep the ventilation openings in the casing of the appliance or in the built-in structure free of obstructions.**
-  **Do not use electrical appliances inside the compartments for storing the frozen food of the appliance, if these are not of the type recommended by the manufacturer.**
-  **If the power cable is damaged, disconnect the appliance from the power supply, contact the technical assistance service for replacement with a model approved by the manufacturer.**

The appliance must only be used within the temperature range specified by the manufacturer. To determine the correct operating range read the **CE plate** in the appliance and locate it on the table given below: Example: model A70/1NE N => 0°+10°C

Series	Temperature
N	0°+10°C
M	-2°+8°C
B	-18°-22°C
T	-10°-30°C
BJ (static/ventilated)	-5°-10°C/-10°-20°C
P (fish)	-6°+4°C

# 2.0 SPECIFICATIONS

<i>made in Italy</i>		Temperature	
Matricola / Serial n°		Modello / Model	
P. Supply	Frequency	Current	
Abs. Power	Clima Class	IP	
Cool. Cap.	Refrigerant	Quantity	
Cool. Cap.	Refrigerant	Quantity	
Cool. Cap.	Refrigerant	Quantity	
2006/95/CE 2004/108/CE		Ins. HFC 245	
			

The appliance conforms to European directives as given in detail on the attached **certificate of EC conformity**. The data are featured on the CE plate inside the appliance.

In airborne noise emissions, the A-weighted sound pressure level is less than 70 dB (A).

Models concerned: A 30-50-60-80-70-120-140 N-M-P-B-T-BJ V

TA-TP-TS-TV-HS 09-13-17-21 M-B-R

## 3.0 INSTALLATION

### 3.1 INTRODUCTION

The appliance must be used solely for the purpose it was specifically designed for, i.e. for storing food within the temperature range stated by the manufacturer. Understand that any other use shall be considered improper.

 **The Manufacturer declines all responsibility in the event of improper use of the appliance.**

### 3.2 DAMAGE CHECK

Check all parts of the appliance for damage and that the standard-issue accessories inside are as expected.

### 3.3 CLEANING

Remove the PVC film covering the inside and outside of the appliance. Clean the compartment inside using a sponge damp with lukewarm water.

### 3.4 WORKING CONDITIONS

Make sure the room in which the appliance is installed meets the following conditions:

- Room temperature in the range 5°C to 43°C.
- Positioning away from sources of heat and in a well ventilated area.
- Leave at least 30mm of space between the back of the cabinet and the wall.
- In case of multiple side-by-side cabinets, leave at least 10mm of space between one cabinet and the other.
- Do not place objects above the cabinet.
- Do not place food on the bottom of the cabinet but place the shelf in the lowest position.

### 3.5 LEVELLING

Move the appliance into its final position and adjust the screw-type feet until you have the unit perfectly level.

### 3.6 PUTTING INTO OPERATION

Before connecting to the power mains, make sure:

- Voltage and frequency are in conformity with the working conditions featured on the CE plate inside the appliance: maximum tolerance is  $\pm 10\%$  of the rated value.
- The power supply circuit meets regulations in force.
- The electrical system is fitted with a residual current circuit breaker (automatic cutout).
- Earthing is compulsory.

If the appliance is going to be switched off for a length of time, it is best to disconnect its power cord.

 **The manufacturer cannot be held responsible for damage or injury resulting from failure to earth the appliance or inefficient earthing thereof, incorrect installation, tampering, poor maintenance and incompetent use, or resulting from failure to comply with electrical safety standards in force in the country where the appliance is used.**

 **For connection to the network, provide a device that ensures complete disconnection in the conditions of overvoltage category III, in accordance with the installation rules.**

## 4.0 CONTROL PANEL

	To see the max. stored temperature; in programming mode it browses the parameter codes or increases the displayed value.		To see the min stored temperature; in programming mode it browses the parameter codes or decreases the displayed value.
	To start a manual defrost.		To display target set point; in programming mode it selects a parameter or confirm an operation.
	Turn on or turn off the light (only versions with lighting)		To switch the instrument off.
+	To lock & unlock the keyboard.	+	To enter in programming mode.
+	To return to the room temperature display.		Fan speed control (Pastry Clima static/ventilated only).

LED	MODE	FUNCTION
	On	Compressor enabled
	Flashing	Anti-short cycle delay enabled
	On	Defrost enabled
	Flashing	Drip time in progress
	On	Fans enabled
	Flashing	Fans delay after defrost in progress.
	On	An alarm is occurring
	On	Continuous cycle is running
	On	Energy saving enabled
	On	Light on
°C/°F	On	Programming phase
	Flashing	Measurement unit

## 5.0 OPERATION

### 5.1 SWITCHING ON

Once you have turned on the appliance's power, the control panel (4.0) automatically runs a lamp test. The LEDs and displays come on for a few seconds, after which the display reads "OFF" (control panel goes into Standby mode).

Make sure the appliance is in standby mode.

Press button .

Display  gives the temperature inside the appliance.

Led  flashes for a few minutes with the compressor off (compressor start delay).

LED  stops flashing and remains steadily lit, the compressor starts and works non-stop to reach the setpoint temperature that has been entered.

At this stage, LED  may switch on and off, depending on the start temperature set for the fans in the general parameters (except for STATIC version).

### 5.2 SWITCHING OFF

 If leaving the appliance switched off for some time, disconnect the power cord.

Press button .

Display  reads "OFF".

Control panel goes into Standby.

### 5.3 SET POINT

The setpoint is the temperature value the appliance is required to reach and maintain, without exceeding it by more than 2°C.

### 5.4 EDITING SET POINT

Press button .

The setpoint value entered flashes on display ; LED  flashes.

Use buttons  and  to increase or decrease the setpoint until you reach the desired.

Press button  to store the value.

### 5.5 WORK CYCLE

Once the temperature setpoint is reached, led  goes off and the compressor stops. It starts again when the temperature rises above said setpoint by 2°C. led  lights again and the cycle is repeated.

Example: if the setpoint is 2°C, the compressor will work to keep the temperature within the range 2°C to 4°C.

## 6.0 SENSORS

---

The appliance features two NTC temperature sensors: one measures the temperature inside the compartment and the other the temperature of the evaporator.

**The second sensor has been fitted so that fan operation is optimized and defrosting times shortened.**

## 7.0 DEFROSTING

---

During operation, frost forms on the surface of the evaporator: this is normal but it reduces evaporator efficiency. The purpose of defrosting is to eliminate frost and restore full efficiency.

### 7.1 AUTOMATIC DEFROSTING

The defrosting process is automatic and takes place at preset intervals.

The purpose of the evaporator sensor is to produce the utmost efficiency and shorten defrosting times by measuring the end-of-defrosting temperature.

The time taken can vary depending on how much frost there is.

### 7.2 MANUAL DEFROSTING

 **This mode should be used only when further defrosting is required in addition to that performed automatically by the appliance.**

Make sure the appliance is working (5.1).

Hold button  down for a few seconds.

LED  comes on.

The defrosting cycle stops automatically. How long it takes depends on how much frost there is on the surface of the evaporator.

Once defrosting has finished, the appliance resumes normal operation.

## 8.0 CONDENSATION EVAPORATION

---

The appliance features a system for the automatic collection and evaporation of water produced by defrosting.

## 9.0 MINIMUM AND MAXIMUM TEMPERATURE

---

You can view the maximum or minimum temperature recorded. This feature allows you to check that products are being stored properly when the business has been closed for a period and hence the appliance has not been monitored.

### 9.1 VIEWING MAXIMUM TEMPERATURE RECORDED

Press key  the display reads "hi" followed by the maximum temperature value recorded.

## 9.2 VIEWING MINIMUM TEMPERATURE RECORDED

Press key  the display reads "Lo" followed by the maximum temperature value recorded.

## 9.3 ERASING RECORDED VALUES

View the maximum or minimum temperature using buttons  or .

Press button .

The "rST" message flashes 3 times to confirm that values have been successfully erased.

# 10.0 OPERATING TIPS

---

When filling with food, if the operation takes more than 5 minutes, switch off the appliance and restart it when you have finished filling.

Only fill freezers with pre-frozen products or with non-frozen products whose mass must not exceed 10% of the total content of the compartment.

Do not use the appliance to cool hot food or food that has just been cooked.

Keep products arranged within the perimeter of the rack so as to allow correct airflow inside the compartment.

The maximum load for each shelf is about 30 kg.

Avoid opening the door frequently when the appliance is defrosting (7.0).

# 11.0 ALARMS AND WARNINGS

---

If an alarm is triggered, the warning is provided by LED , which remains on display  until the alarm condition is over. All alarm messages flash, alternating with the temperature, except for P1, which flashes constantly.

## 11.1 LIST OF ALARMS

"P1" flashing; **compartment sensor error.**

- Tells you that compartment sensor is faulty or is not connected properly.
- Emergency mode is triggered automatically, which handles the compressor start and stop cycles based on values set in the general parameters.
- If the alarm persists, contact technical assistance department .

"P2" flashing; **evaporator sensor error.**

- Tells you that evaporator sensor is faulty or is not connected properly.
- If the alarm persists, contact technical assistance department.

"HA" flashing; **high temperature.**

- Tells you that temperature inside the compartment is higher than maximum permissible value.
- The alarm ceases automatically when temperature falls back inside the normal range, when defrosting starts or when the door is opened.
- If the alarm persists, contact technical assistance department.

"LA" flashing; **low temperature.**

- Tells you that temperature inside the compartment is lower than minimum permissible value.
- The alarm ceases automatically when temperature falls back inside the normal range, when defrosting starts or when the door is opened.
- If the alarm persists, contact technical assistance department.

"EE" flashing; **memory anomaly.**

- Anomaly in control panel data. Switch the appliance off and back on again to reset the instrument.
- If alarm persists, contact assistance department.

"dA" flashing; **door open.**

- Tells you that the door has been open too long. Ceases automatically when door is closed.

You can silence the buzzer announcing an alarm by pressing any button.

## 12.0 CONNECTION DIAGRAM

---

- Appendix A
- Appendix B for cabinets Fast Master
- Appendix C for remote cabinets

## 13.0 MAINTENANCE AND CLEANING

---

Routine maintenance work can be carried out by non-specialized personnel, following the instructions given below to the letter.

-  **Before performing any maintenance or cleaning work, disconnect the appliance from the power mains.**

### 13.1 CLEANING THE APPLIANCE

You can wash the inside and outside using a sponge damp with lukewarm water and neutral detergents without chlorine content and without any other chemically aggressive component. Once you have finished cleaning, dry with a soft, dry cloth. You are advised to apply polishing products on outer surfaces only.

-  **Do not wash the appliance with jets of water. Do not use scouring pastes or steel wool. Do not spray detergents on internal and external surfaces.**

### 13.2 CLEANING THE CONDENSER

To keep the appliance working efficiently at all times, clean the condenser at regular intervals. To do this, use a brush with soft bristles or a vacuum cleaner, being careful not to bend the aluminium fins.

-  **The condenser features sharp edges. Wear protective gloves when cleaning.**

### 13.3 CLEANING THE GASKETS

To prevent mold formation, it is recommended the periodic cleaning of the gaskets. Use only lukewarm water or a solution of water and baking soda. Do not use vinegar, detergents, degreasers, or acid products.

-  **Do not use vinegar, detergents, degreasers, or acid products.**

# 14.0 GAS SAFETY DATA SHEET

## 14.1 R290

GWP : 3  
ODP : 0

- Hazards identification: Liquefied gas - Extremely flammable
- First aid measures:
  - Inhalation : Move victim to uncontaminated area by wearing self-contained breathing apparatus. Keep the patient relaxed and warm. Call a doctor. Proceed with artificial respiration if breathing stops.
  - Contact with the skin : In the case of frostbite, spray with water for at least 15 minutes. Apply a sterile gauze. Get medical assistance.
  - Eye contact : Immediately flush eyes with water for at least 15 minutes.
  - Ingestions : Ingestion is considered an unlikely route of exposure..

## 14.2 R452A

GWP : 2141  
ODP : 0

- Composition :

○ Difluoromethane	(R32)	21,9%
○ 2,3,3,3-Tetrafluoropropene	(R-1234yf)	27,2%
○ Pentafluoroethane	(HFC 125)	50,9%
- First aid measures: :
  - Inhalation : in In high concentration it can cause asphyxia. Symptoms may include loss of mobility and / or knowledge. Victims may not realize the asphyxia. Move victims to a ventilated area and keep them warm wearing the self-contained breathing apparatus. Call a doctor. Only perform artificial respiration if the breath has stopped.
  - Contact with skin : Contact with the evaporating liquid may cause skin freezing.
  - Contact with eye : Immediately flush eyes with water. Remove any contact lenses if it is easy to do. Continue to rinse. Rinse with plenty of water for at least 15 minutes. Immediately call for a medical examination. If medical assistance is not immediately available, rinse for another 15 minutes.
  - Ingestion : Ingestion is considered an unlikely route of exposure.

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# 2.0 TECHNISCHE DATEN

Matricola / Serial n°		Modello / Model		Temperature	
P. Supply		Frequency		Current	
Abs. Power		Clima Class		IP	
Cool. Cap.		Refrigerant		Quantity	
Cool. Cap.		Refrigerant		Quantity	
Cool. Cap.		Refrigerant		Quantity	
2006/95/CE 2004/108/CE			Ins. HFC 245		
					

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